



IBC Policy: Testing Lentiviral Vectors for Replication Competent Virus

The lentiviral vector (HIV) preps used in animals or for experiments involving sharps must be evaluated for the appearance of replication-competent retrovirus (RCR) using the three independent assays described below:

1. **Tat-transfer assay** – This assay is based on a reporter HeLa-CD4-LTR- β gal cells containing one integrated copy of the HIV-1 LTR (nts -138 through +80) linked to the β -galactosidase gene. The reporter cell line is highly susceptible to infection. In the case of viral genome recombination that results in the reconstitution of replication-competent HIV-1, the recombined vector will be capable of generating functional tat protein. The tat-expression will lead to activation of the viral LTR-promoter driven the expression of β -galactosidase gene of the reporter cell line. The assay sensitivity has been determined to be as low as 20 tat-transducing units per ml of test medium. The assay is performed as follows: The cells transduced with lentiviral vector are serially passaged three-four times (about 2 weeks), after which the supernatant of the cells is transferred to a reporter HeLa-CD4-LTR- β gal cells. The HIV-1-tat activity is determined by X-Gal (5-bromo-4-chloro-3-indolyl-b-D-galactopyranoside) staining. By this method, vector preparation is considered helper negative when no expression of β gal is detected.
2. **Gag-transfer assay** – This assay is based on the detection of p²⁴gag-protein of the virus in conditioned media obtained from vector-transduced cells. The cells are serially passaged three-four times (about 2 weeks), after which the supernatant of the cells is collected for assessing the level of p24gag by ELISA (p24 ELISA kit, NIH). The detection limit of this method is ≥ 100 pg/mL of p²⁴, which is about 10^3 copies of viral genome per mL. By this method, vector preparation is considered helper negative when p24 concentrations are below detection levels.
3. **Marker-rescue assay** – This assay is based on the direct detection of GFP or other reporter following the transfer. Viral vector stocks are assessed as follows. The cells transduced with lentiviral vector harboring a reporter (GFP) are serially passaged three-four times (about 2 weeks), after which the supernatant of the cells is collected and transferred to HEK293T cells cultured in a 10-cm plate. Seventy-two hours post-transduction, the cells are scored for a reporter expression. Vector stock was considered helper free when no reporter is detected.

References:

1. Charneau, P., G. Mirambeau, P. Roux, S. Paulous, H. Buc, and F. Clavel. 1994. HIV-1 reverse transcription. A termination step at the center of the genome. *J. Mol. Biol.* 241:651–662.
2. Kimpton J, Emerman M. Detection of replication competent and pseudotyped HIV with a sensitive cell line based on activation of an integrated beta-galactosidase gene. *J Virol* 66:2232-2239, 1992.
3. Naldini, L., U. Blomer, P. Gallay, D. Ory, R. Mulligan, F. H. Gage, I. M. Verma, and D. Trono. 1996. In vivo gene delivery and stable transduction of nondividing cells by a lentiviral vector. *Science* 272:263–267.